











በጣም በርካታና የተለያዩ የህክምና ምርመራዎች እሉ፤ በመሆኑም የተለያዩ አዘገጃጀቶችን ይጠይቃሉ። ለምሳሌ፣ አንዳንዶቹ ምርመራዎች መጸምን የሚጠይቁ ሲሆን አንዳንዶቹ ግን አያስፈልጋቸውም። ለእርስዎ ስለሚሰጣቸው የምርመራ ዓይነት ይክተርዎን ያነጋግሩ።

የሰገራ ናሙና ኣወሳሰድ

<p>ነጩን ወረቀት ይዘርጉት። ከሽንት ቤቱ የመጸዳጃ ገንዳ ውስጥ በውሀው ላይ ያስቀምጡት።</p>	
<p>ሰገራ ይውጡ። የሚወጣው ሰገራ በወረቀቱ ላይ ያርፋል። አረንጓዴውን ክዳን ጠምዘዝ አድርገው ይክፈቱት። እንጨቷን ያውጧት።</p>	
<p>ወረቀቱ ወድቆ ሰገራው ውሀ ከመንካቱ በፊት የሰገራውን ናሙና ጠንቆል አድርገው ያንሱ። ሠገራውን በትንሿ እንጨት ነገር ዛቅ አድርገው ያንሱ።</p>	
<p>የተባረባረው የእንጨቷ ጫፍ በሰገራው ናሙና መሞላቱን ያረጋግጡ።</p>	
<p>እንጨቷን መልሰው ብልቃጧ ውስጥ ያስገቧት። አረንጓዴዎን የብልቃጥ ክዳን በደንብ ጥምዘዝ አድርገው ይዝንት። እንደገና አይክፈቱ። በመጸዳጃው ገንዳ ውሀ ይልቀቁበት። ነጩ ወረቀት የቧንቧ ስራው ላይ ችግር አይፈጥርም።</p>	
<p>ብልቃጧንና የላቦራቶሪ ወረቀቱን በፕላስቲክ ቀረጫቱ ውስጥ አድርገው ወደላቦራቶሪው መልሰው ይላኩት።</p>	

There are many different kinds of medical tests and they require different kinds of preparation. For example, some tests require fasting while other tests do not. Talk to your doctor about the test preparation that is right for you.

Stool Sample Collection

<p>Unfold white paper. Put inside toilet bowl on top of the water.</p>	
<p>Have a bowel movement. Stool will be on top of the paper. Twist and lift off the green cap. Pull out the stick.</p>	
<p>Take stool sample before paper falls and stool touches water. Scrape stool with the stick.</p>	
<p>Make sure the grooved part of the stick is covered with stool sample.</p>	
<p>Put stick back in bottle. Screw the green cap on tight. Do not reopen. Flush toilet. White paper will not harm plumbing.</p>	
<p>Put bottle and lab requisition in plastic bag and return to lab.</p>	